



metacarpal wrist guard

knee guard

wrist guard with splint

elbow guard for epicondylitis

soft cervical collar

bandage maximum aeration

chest belt

elbow guard

shoe for plaster

meta thumb wrist guard

ankle guard

elastic band for ankle

knee guard

arm support bag

elastic compression bandages

bulvar wool bandage

short tight with strong compression

velcro seal bandage

strong compression tight

Intex

CATALOGUE**PRODUCTS**

BEIGE RANGE	4
MALLAFIX	6
MALLAFIX PLUS	8
FORTIFLEX	9
CHEST BELT RANGE	10
BLACK RANGE	11
BLUE RANGE	12
WRIST AND ELBOW GUARDS	14
PLASTER	15
WOOL/COTTON BANDAGES	16
THERMAL ORTHOSIS	17
RENGO	18
BENTERM	19
STOCKINGS RANGE	20
OTHERS RANGE	22

This style is particularly indicated for injuries of the joints, sprains, chronic ailment, traumatism, etc and as a preventive measure against efforts or long exposures to vibration and overloads both in sports and occupational.

- Bi elastic orthosis in anatomical fabric.
- Variable mesh thickness.
- Anatomic shapes.

7.133 / Knee Guard *with lateral flexible reinforcement.*

SIZE	MEASUREMENT (cm.)
Small	30,0 to 33,5
Medium	33,5 to 37,0
Large	37,0 to 40,5
Extra	40,5 to 44,0

Measures indicate the perimeter of the knee without bending.



ART. 7.133

7.135 / Metacarpal wrist guard *with splint.*

SIZE	MEASUREMENT (cm.)
Small	13,0 to 16,0
Medium	17,0 to 20,0
Large	21,0 to 24,0

Measures indicate the perimeter of the wrist.

Valid for right or left hand.

With reinforcing tape for extra support and immobilization.



ART. 7.135



7.137 / Elastic band *for wrist.*

ONE SIZE FITS ALL



Elastic bands *for wrist and ankle.*

- Elastic bands self attachable at any point of themselves.
- They allow to configure different kinds of compressive dressings.
- Fastening clips not required.
- They can be used individually or as a reinforcement of an orthosis.

7.138 / Elastic band *for ankle.*

ONE SIZE FITS ALL



7.134 / Ankle guard

SIZE	MEASUREMENT (cm.)
Small	18,0 to 20,5
Medium	20,5 to 23,0
Large	23,0 to 25,5
Extra	25,5 to 28,0

Measures indicate the perimeter of the ankle above the malleolus.



7.136 / Enclosed wrist guard

SIZE	MEASUREMENT (cm.)
Extra small	12,0 to 14,0
Small	14,0 to 16,0
Medium	16,0 to 18,0
Large	18,0 to 20,0
Extra	20,0 to 22,0

Measures indicates the perimeter of the wrist.

Mallafix orthosis elastic supports have been designed to protect the joints.

The absence of rigid splints allows to maintain a controlled mobility.

Especially recommended for joint injuries and as a preventive measure to prolonged efforts, vibration and overloads.

Based in multimesh technology and bi elastic fabric.

Great adaptability and ease of use.



7.097 / Enclosed metacarpal wrist guard

- Protects carpal and metacarpal allowing mobility.
- Available in white and blue colour.



SIZE	MEASUREMENT (cm.)
Extra Small	12,0 to 14,0
Small	14,0 to 16,0
Medium	16,0 to 18,0
Large	17,5 to 20,0
Extra Large	19,0 to 21,5

Measures indicate wrist perimeter.



ART. 7.097

7.098 / Elbow guard


- High efficacy against lateral and medial epicondylitis.
- Available in white colour.

SIZE	MEASUREMENT (cm.)
Small	23,0 to 26,5
Medium	25,0 to 29,0
Large	27,5 to 31,0
Extra Large	29,5 to 33,5

Measures indicate the elbow perimeter without bending.



ART. 7.098



7.094 / Lumbar bandage

- Front height: 14 cms.
- Back height: 24 cms.
- Available in white colour.

SIZE	MEASUREMENT (cm.)
Small	68,5 to 84,0
Medium	84,0 to 100,0
Large	99,0 to 115,5

Measures indicates waist perimeter.



7.095 / Ankle guard

- Designed to avoid wrinkles in the instep area.
- Available in white and blue colour.

SIZE	MEASUREMENT (cm.)
Small	19,0 to 21,5
Medium	21,0 to 24,0
Large	23,0 to 26,5
Extra Large	25,0 to 29,0

Measures indicates the perimeter of the ankle above the malleolus.



7.096 / Knee guard

- Designed to avoid wrinkles behind the knee (plopital area).
- Available in white and blue colour.

SIZE	MEASUREMENT (cm.)
Small	29,5 to 33,5
Medium	31,5 to 36,0
Large	33,5 to 38,5
Extra Large	35,5 to 41,0

Measures indicates the knee perimeter without bending.

Bi elastic orthosis similar to Mallafix range products with silicone cushion for a bigger therapeutic effect.

7.120 / Lumbar bandage

- With silicone cushion in the lumbar zone at the sacrum level.

SIZE	MEASUREMENT (cm.)
Small	68,5 to 84,0
Medium	84,0 to 100,0
Large	99,0 to 115,5

Measures indicate the perimeter of the waist.



ART. 7.120

7.121 / Ankle guard

- With two silicone cushions in the retromalleolar hollows.

SIZE	MEASUREMENT (cm.)
Small	19,0 to 21,5
Medium	21,0 to 24,0
Large	23,0 to 26,5
Extra Large	25,0 to 29,0

Measures indicates the ankle perimeter above the malleolus.



ART. 7.121

7.122 / Knee guard

- With two lateral reinforcement.
- Incorporates a silicone ring around the kneecap.

SIZE	MEASUREMENT (cm.)
Small	29,5 to 33,5
Medium	31,5 to 36,0
Large	33,5 to 38,5
Extra Large	35,5 to 41,0

Measures indicate the knee perimeter without bending.



ART. 7.122


FORTIFLEX

7.123 / Sacro lumbar bandage

- Front height between 17 and 19 cms.
- Back height 26,5 cms.
- Reinforcement bands at the back for a higher compression in the lumbar area.
- With a handle for easy closing.
- Incorporates 6 self adaptable flexible splints.

SIZE	MEASUREMENT (cm.)
Small	69 to 84
Medium	80 to 95
Large	91 to 106
Extra Large	102 to 117

Measures indicate the waist perimeter.



ART. 7.123

CHEST BELT RANGE

Intex

6.160 / One band chest belt (16 cms)

- Breathable fabric.

SIZE	MEASUREMENT (cm.)
Small	65 to 90
Medium	75 to 100
Large	85 to 110



2, 3 and 4 band chest belts

- Chest belts made with 8 cms bands linked together with different heights: 16, 24, 32 cms.
- The way the bands are linked together allows adaptability without producing wrinkles.



7.077 / 7.070 Double band chest belt

SIZE

ONE SIZE FITS ALL (7.077) for man
ONE SIZE FITS ALL (7.070) for woman

7.077 3B / 4B 3 and 4 band chest belt

SIZE	MEASUREMENT (cm.)
1	75 to 95
2	85 to 105
3	100 to 120
4	110 to 130



7.072 / Bandage maximum aeration

- Lumbar breathable bandage. Height: 27 cms.

SIZE	MEASUREMENT (cm.)
1	62 to 82
2	75 to 95
3	88 to 107
4	100 to 120



7.128 / Black lumbar bandage

- Black lumbar bandage with Velcro seal and lateral adjustable splints allowing a customized compression adjustment depending on the activity at any time.
- With 4 self adaptable flexible plastic splints in the lumbar zone.
- Front height: 17 cms.
- Back height: 27 cms.
- Made with an elastic material that offers a soft touch and an easy anatomical adaptation.
- Available in 3 sizes.

SIZE

Small	for waist perimeters between 80 and 94 cms
Medium	for waist perimeters between 95 and 109 cms
Large	for waist perimeters between 110 and 125 cms



Punctual and repetitive efforts



Prolonged periods of driving



Sports activities



Prolonged exposure to knocks or vibrations



Weightlifting



ART. 7.128

The blue is a whole range of orthosis products for the upper limb, made with anatomical bi elastic knitting.

7.126 / Short thumb abductor wrist

- To keep the thumb in abduction for injuries in the metacarpophalanx joint of the thumb.
- With a removable lateral splint and manually malleable.

SIZE	MEASUREMENT (cm.)
Small	13 to 16
Medium	17 to 20
Large	21 to 24

Measures indicates the perimeter of the wrist.



ART. 7.126

7.125 / Meta-thumb wrist guard

- Metacarpal wrist guard with splint and thumb abductor.
- High degree of stabilization and immobilization.
- Similar to 7124 but with a thumb abductor and a second splint.

SIZE	MEASUREMENT (cm.)
Right - Small	13 to 16
Right - Medium	17 to 20
Right - Large	21 to 24
Left - Small	13 to 16
Left - Medium	17 to 20
Left - Large	21 to 24

Measures indicates the perimeter of the wrist.



ART. 7.125



7.127 / Elbow guard for epicondylitis

- Recommended for lateral and medial epicondylitis.
- Incorporates two removable metal splints allowing for precise positioning.

SIZE	MEASUREMENT (cm.)
Small	24 to 30
Large	27 to 34



7.124 / Metacarpal wrist guard with splint

- Protection and stabilization of the bottom radio ulnar joint of the carpus and metacarpus.
- High degree of stabilization and immobilization.
- Manually malleable and removable palmar splint.

SIZE	MEASUREMENT (cm.)
Right - Small	13 to 16
Right - Medium	17 to 20
Right - Large	21 to 24
Left - Small	13 to 16
Left - Medium	17 to 20
Left - Large	21 to 24

Measures indicates the wrist perimeter.

7.006 / Wrist guard. *Adjustable lock.*

- Width: 8 cms.
- Incorporates an elastic ring to facilitate its placement.

ONE SIZE FITS ALL



ART. 7.006

7.035 / Anti-epicondylitis bandage

- Recommended for lateral and medial epicondylitis.
- Incorporates two removable metal splints allowing for precise positioning.

SIZE	MEASUREMENT (cm.)
Small	24 to 30
Large	27 to 34



ART. 7.035

7.104 / Wrist guard. *Metacarpian with splint.*

- Protection and immobilization of the bottom radio ulnar joint of the carpus and metacarpus.
- High degree of stabilization and immobilization.
- Manually malleable and removable palmar splint.

SIZE	MEASUREMENT (cm.)
Right - Small	13 to 16
Right - Medium	17 to 20
Right - Large	21 to 24
Left - Small	13 to 16
Left - Medium	17 to 20
Left - Large	21 to 24



ART. 7.104



7.048 / Shoe for plaster

- Anti-slip and waterproof sole.
 - Applies to both left and right foot.
 - For plaster foot or post operative processes.
-

SIZE

Small
Medium
Large



7.081 / Arm support bag

- Recommended to those who have to carry the arm in plaster.
 - Provides support to the arm and prevents the hand is falling.
-

ONE SIZE FITS ALL



7.080 / Arm support sling

- Recommended to those who have to carry the arm in plaster.
 - Provides a buckle for adjusting the height of the arm.
-

SIZE

Adult
Children

WOOL/COTTON BANDAGES

Intex

Traditional bandages that provide an effective protection for the waist. Provide a slight compression and a constant temperature. High degree softness.

7.017 / Wool bandage. Adjustable lock.

- Size 1: 27 cms height and size: 33 cms height.

SIZE	MEASUREMENT (cm.)
1	60 to 90
2	90 to 120

7.087 / Tubular wool bandage

- Available in 27 and 33 cms height.

SIZE	MEASUREMENT (cm.)
27/5	60 to 75
27/6	75 to 90
27/7	90 to 105
27/8	105 to 120
33/5	60 to 75
33/6	75 to 90
33/7	90 to 105
33/8	105 to 120

7.086 / Velcro seal cotton bandage

- Size 1: 26 cms height; sizes 2 and 3: 30 cms height.

SIZE	MEASUREMENT (cm.)
1	70 to 88
2	88 to 114
3	110 to 130

7.086 / Tubular cotton bandage

- Small and medium size: 26 height.
- Large and extra large: 30 cms height.

SIZE	MEASUREMENT (cm.)
Small	75 to 82
Medium	82 to 88
Large	88 to 95
Extra Large	95 to 105



Superfine 100% merino wool styles and elastic yarns. Provides a nice thermal protection. High capacity to absorb moisture. Light compression. Excellent thermal regulator. Do not cause itching.

7.150 / Wool bandaje merino

SIZE	MEASUREMENT (cm.)	HEIGHT (cm.)
Small	60 to 75	27
Medium	75 to 90	27
Large	90 to 105	33
Extra	105 to 120	33

Measures indicate the waist perimeter.



7.151 / Knee guard merino

SIZE	MEASUREMENT (cm.)
Small	29,5 to 33,5
Medium	31,5 to 36,0
Large	33,5 to 38,5
Extra	35,5 to 41,0

Measures indicates the knee perimeter without bending.



7.152 / Elbow guard merino

SIZE	MEASUREMENT (cm.)
Small	23,0 to 26,5
Medium	25,0 to 29,0
Large	27,5 to 31,0
Extra	29,5 to 33,5

Measures indicate the elbow perimeter without bending.



Made of cotton thread, with an anatomical shape to offer uniform compression for the treatment of different pains which require compression.

PRODUCT	SIZE	MEASUREMENT (cm.)
7.089 / ENCLOSED WRIST GUARD	1	16 to 18
	2	17 to 19
	3	18 to 20
	4	19 to 21
	5	20 to 22
7.090 / KNEE GUARD	2	25
	3	26
	4	27
	5	28
	6	29
	7	30
	8	31
	9	32
	10	33
	11	34
	12	35
	13	36
	14	37
7.091 / ANKLE GUARD	2	13 to 17
	3	14 to 18
	4	15 to 19
	5	16 to 20
	6	17 to 21
	7	18 to 22
	8	19 to 23
	9	20 to 24
	10	21 to 25
	11	22 to 26
	12	23 to 27
	13	24 to 28
	14	25 to 29
7.092 / SHORT TIGHT WITH STRONG COMPRESSION	2	30
	3	31
	4	32
	5	33
	6	34
	7	35
	8	36
	9	37
	10	38
	11	39
	12	40
	13	41
	14	42
7.093 / METACARPAL WRIST GUARD	1	14,5 to 15,5
	2	15,5 to 16,5
	3	16,5 to 18,0
	4	18,0 to 19,5
	5	19,5 to 21,0

Enclosed wrist guard



ART. 7.089

Knee guard




ART. 7.090

Ankle guard



ART. 7.091

Short tight with strong compression



ART. 7.092

Metacarpal wrist guard



ART. 7.093

Made with thermogenic fibers with a high calorific value. Soft elasticity and pleasant touch. Wide range of colours.



ART. 7.088

7.088 / Wrist guard. Adjustable velcro seal.

- Elastic support with a velcro seal.
- Width: 6 cms.
- Elastic ring for a correct positioning.

SIZE

ONE SIZE FITS ALL - Beige
ONE SIZE FITS ALL - White



ART. 2.761

2.761 / Bandage. Adjustable velcro seal.

- Every size is 30 cms height.
- Available in white colour.

SIZE

MEASUREMENT (cm.)

Small	75 to 93
Large	87 to 120
Super	108 to 150



ART. 2.760

2.760 / Tubular Bandage.

SIZE


MEASUREMENT (cm.)


HEIGHT

Beige, 1	60 to 72	26
Beige, 2	73 to 85	28
Beige, 3	86 to 98	30
Beige, 4	99 to 111	32
Beige, 5	112 to 120	34
Beige, 6	121 a 140	36

White, 1	60 to 72	26
White, 2	73 to 85	28
White, 3	86 to 98	30
White, 4	99 to 111	32
White, 5	112 to 120	34
White, 6	121 to 140	36

Blue, 1	60 to 72	26
Blue, 2	73 to 85	28
Blue, 3	86 to 98	3
Blue, 4	99 to 111	32
Blue, 5	112 to 120	34
Blue, 6	121 to 140	36





ELASTIC STOCKINGS WITH DECREASING COMPRESSION.

Tights with maximum compression in the ankle and minimum compression in the thigh (decreasing) for the prevention or treatment of venous insufficiency.

They attenuate the risk factors and curb the evolution of the illness.

We have different levels of compression, light compression (40 Den), normal compression (140 Den) and strong compression (280 Den).

Compression tights should be put on as soon as we get up, as later on legs begin to swell, making it difficult to put them on. This is specially important in strong compression.

RECOMMENDATIONS TO CHOOSE COMPRESSION:

COMPRESSION	DEN.	PRESSURE	INDICATIONS
LIGHT	40 - 70	15-21 mm Hg.	Long trips.
NORMAL	140	22-29 mm Hg.	Pregnancies. Risk of venous insufficiency in jobs that require standing.
STRONG	280	30-40 mm Hg.	Prevention of varicose veins. Mild chronic venous insufficiency. Trombophlebitis. Symptomatic varicose veins.

Measurement table (cm.)

Size	Ankle	Calf	Thigh
2	18/19	30/31	41/42
3	20/22	32/34	43/45
4	23/25	35/38	46/49
5	26/28	39/41	50/53
6	26/28	42/44	54/59

ART.	NAME	SIZE	COLOUR
7.203	STRONG COMPRESSION TIGHT	2	Beige
		3	Beige
		4	Beige
		5	Beige
		6	Beige
7.203	STRONG COMPRESSION TIGHT, PREGNANCY	2	Beige
		3	Beige
		4	Beige
		5	Beige
7.203	LONG STOCKING, STRONG COMPRESSION, LACE	2	Beige
		3	Beige
		4	Beige
		5	Beige
7.203	STRONG COMPRESSION SHORT STOCKING	2	Beige
		3	Beige
		4	Beige
		5	Beige
7.022	NORMAL COMPRESSION TIGHT	2	Beige
		3	Beige
		4	Beige
		5	Beige
		6	Beige
7.022	NORMAL COMPRESSION SHORT STOCKING	2	Beige
		3	Beige
		4	Beige
		5	Beige
7.201	SHORT STOCKING 40 DEN.	ONE SIZE FITS ALL	Black
		ONE SIZE FITS ALL	Beige
		ONE SIZE FITS ALL	Brown
		ONE SIZE FITS ALL	Grey
7.201	MEN SOCK 40 DEN.	ONE SIZE FITS ALL	Black

7.079 / Elastic compression bangades

- High compression elastic bandages for the treatment of venous insufficiency or for the compressive bandage of the joints.
- Clips for attachment provided.
- The length of the 3 sizes is 2 meters resting and 5 meters extended.
- The width of each size are 5, 7 and 10 cms.

SIZE

5 x 5

7 x 5

10 x 5



ART. 7.079

7.550 / Soft cervical collar

- Foam collar for using in cases of cervical whiplash, dislocations or bruises on the neck area.
- Its use keeps the neck and the cervical spine straight.

SIZE

PATIENT'S MEASUREMENT (cm.)

Small	36 to 41
Medium	41 to 48
Large	46 to 51



ART. 7.550

NOTES

Intex, s.a.

Paseo Lizarkola s/n.
Apdo. Correos 22
20140 Andoain (Guipúzcoa)
España

Tel.: 00 34 943 590 900 / 00 34 943 590 851
Fax: 00 34 943 591 400

intex@intex-sa.com

www.intex-sa.com

Location: lat. 43° 12' 43,2 N
lon. 02° 00' 55,9 W